MQ XL Carbon Stroller www.mobiquip.com

MQ
MobiQuip

Table of Contents

1. Features
2. Warning
3. Measurement
4. Assemble
5. Function

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A circular icon containing the letters 'EN' in a bold, sans-serif font.

WARNING!

READ AND FOLLOW ALL WARNINGS AND INSTRUCTIONS! FAILURE TO DO SO COULD RESULT IN SERIOUS INJURY OR DEATH!

NEVER leave child unattended. ALWAYS keep child in view while in the stroller.

Avoid serious injury from falling or sliding out.

ALWAYS use harness.

When making adjustments to the stroller, ALWAYS ensure that all parts of the child's body are clear of any moving parts of the stroller.

DO NOT use the stroller as a seat for the occupant when in a motor vehicle.

DO NOT use stroller if any parts are missing or broken.

DO NOT allow child to climb on or around stroller.

DO NOT allow child to stand on the footrest. stroller may tip over and injure child.

NEVER leave stroller on a hill or incline, even with brakes set, stroller may tip injuring the child.

To avoid accidental rolling, ALWAYS set brakes when stroller is not being pushed.

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EN

Measurement	
Width	55cm
Height	112cm
Length	96cm
Footrest height	112cm
Footrest dimensions	44 x 19cm
Seat width	40cm
Seat depth	27cm
Back height	56cm
Back width	40cm
Product weight	7.1kg
Weight capacity	60kg

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Plegar el cochecito

EN**Product Identification**

Your stroller is equipped with a frame label, which is located on the frame (under the seat) of your stroller. The frame label contains information about your stroller, for example the product name, manufacturer name, the name of the authorised representative(s), Unique Device Identifier (UDI) and technical data of your stroller. Below are example frame labels, the information on the frame label of your stroller will correspond to your specific stroller. Van Os Medical UK Limited is the importer of this product for the UK. This product is manufactured by Lerado Zhongshan Peaceful Cove Business Trading Co., Ltd, who are represented in the United Kingdom by UK APLAN CORPORATION LTD. Mobiquip Limited is the exclusive distributor of this mobility product. The manufacturer's product code for your stroller is MK1008.

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MQ XL Carbon Stroller
Paediatric travel buggy for indoor and outdoor use

Lerado Zhongshan Peaceful Cove
Business Trading Co., Ltd.
No. 26 Guangfu Road, Xiaolan Town,
Zhongshan City, Guangdong Province,
528414, China

IMPORTER FOR EU27:

Van Os Medical B.V.
Koperslagerij 3, 4651 SK Steenbergen
The Netherlands

IMPORTER FOR UK:

Excel Mobility Ltd.
Excel House, Ashbrooke Park,
Sherburn-In-Elmet, North Yorkshire,
LS25 6PJ, United Kingdom

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011680001

Product Identification

Van Os Medical B.V. is the importer of this product for the EU27. This product is manufactured by Lerado Zhongshan Peaceful Cove Business Trading Co., Ltd, who are represented in the European Union by Y. SUNGHANDELSVERTRETUNG. Mobiquip Limited is the exclusive distributor of this mobility product. The manufacturer's product code for your stroller is MK1008.

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MQ XL Carbon Stroller

Paediatric travel buggy for indoor and outdoor use

Lerado Zhongshan Peaceful Cove
Business Trading Co., Ltd.
No. 26 Guangfu Road, Xiaolan Town,
Zhongshan City, Guangdong Province,
528414, China

IMPORTER FOR EU27:

Van Os Medical B.V.
Koperslagerij 3, 4651 SK Steenberghe
The Netherlands

IMPORTER FOR UK:

Excel Mobility Ltd.
Excel House, Ashbrooke Park,
Sherburn-In-Elmet, North Yorkshire,
LS25 6PJ, United Kingdom

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+

011680001

MobiQuip Ltd
The Tower, Pera Business Park,
Nottingham Rd,
Melton Mowbray
LE13 0PB

Lerado Zhongshan Peaceful Cove
Business Trading Co., Ltd.
No. 26 Guangfu Road, Xiaolan Town,
Zhongshan City, Guangdong Province,
528414, China

IMPORTER FOR UK: (For UK customers)
Excel Mobility Limited
Excel House, Ashbrooke Park, Sherburn-
In-Elmet, North Yorkshire, LS25 6PJ,
United Kingdom

IMPORTER FOR EU27: (For EU27 customers)
Van Os Medical B.V.
Koperslagerij 3,
4651 SK, Steenberghe,
The Netherlands

+ 44 (0) 1977 681 400

customerservice@excelmobility.uk

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UK Aplan Corporation Ltd
KOMAS Medical Technology GmbH